GyneFIX[®]

Physician's Information
Information pour le médecin
Informatie voor de arts
Information für den Arzt
Informazioni per il medico
Información para médicos

CE 0120

MANUFACTURER - FABRICANT - FABRIKANT
HERSTELLER - PRODUZIONE - FABRIKANTE:

*Controlled Release
for the Enhancement
of Quality of Life
of Women*

CONTREL EUROPE NV

Incubation & Innovation Center,
Technology Park, University Ghent,
9052 Zwijnaarde (Ghent), Belgium.
info@control.be

Contact: +32 476 885 738
only for insertion related questions

I. Name of Product

GyneFix® Intrauterine Copper Contraceptive for interval and postabortal use (< 10 weeks).

II. Composition

Qualitative composition: Cuprum 99.99% (high purity copper)
Quantitative composition: 420 mg (standard version); 280 mg (mini version)

III. Description

GyneFix® is a new concept in intrauterine device (IUD) technology. It is fixed to the uterine fundus, is frameless and entirely flexible. GyneFix® consists of 6 copper sleeves (standard version) or 4 copper sleeves (mini version), each 5 mm long and 2.2 mm in diameter, threaded on a length of polypropylene suture material. The sleeves are prevented from sliding off the material by the upper and lower sleeves, being crimped onto the thread. The proximal end of the thread is provided with a knot which at insertion is placed in the fundal myometrium with an inserter for anchoring the device. The total surface of copper, including the inner and outer surfaces, is 330 mm² for the standard version and 200 mm² for the mini version. This new implant device has no plastic body, making it a completely flexible unit. GyneFix® is preloaded in a sterile inserter.

IV. Clinical Pharmacology

Available data indicate that the contraceptive effectiveness of copper IUDs is enhanced by a minute quantity of copper being released continuously into the uterine cavity. The exact mechanism by which copper enhances the contraceptive effect of an IUD has not been conclusively demonstrated. Various hypotheses have been advanced; the major one being that copper placed in the uterus interferes with enzymatic or other processes that regulate blastocyst implantation. In addition, animal studies suggest that copper may play a role in reducing sperm transport within the uterine environment.

V. Indication and Clinical Use

Indication

GyneFix® is indicated for contraception.

Clinical Use

A candidate for use of GyneFix® is any healthy, non-pregnant, female.

Time of insertion

Experience has shown that GyneFix® can be inserted at any time during the menstrual cycle. It is evident, however, that the physician should rule out the presence of a pregnancy if insertion of the device is planned in the second half of the menstrual cycle.

Clinical trials have shown that GyneFix® can be inserted immediately after spontaneous or induced first trimester abortion of less than 10 weeks gestational age. GyneFix® is unsuitable for use in the immediate postpartum period. A period of at least 6 weeks should be allowed until the uterus has regained its normal size or the subject has had at least one normal menstrual period. In breastfeeding women, however, a period of 8 to 12 weeks following full term delivery is recommended because of the higher risk of perforation reported in some studies with traditional IUD's.

VI. Use-Effectiveness Studies

GyneFix® has been studied since 1985 and over 15,000 women years of experience (international, multicenter, non-comparative and comparative clinical trials, including a large proportion of nulligravida/nulliparous women) have been collected.

1. Effectiveness of GyneFix®

GyneFix® is a very effective contraceptive device. The pregnancy rates are comparable or lower when compared with the most effective currently used high-load copper devices, and to oral contraceptives and sterilization.

Use in young parous and nulligravida/nulliparous women

Studies have shown that young parous and nulligravida/nulliparous women using GyneFix® are no more susceptible to pregnancy than other age groups, the annual pregnancy rates being close to zero.

The low pregnancy rate with GyneFix® is attributed to the high delivery of the copper ions in the uterine cavity. The mini version is more suitable for small uterine cavities.

2. Bleeding/pain with GyneFix®

Due to its design characteristics, flexibility, absence of a frame and the reduced space it occupies within the uterine cavity, GyneFix® is highly tolerated and the removal rate for pain is very low.

The removal rate for bleeding is also low with GyneFix® and is ascribed to the small surface area and the characteristics of the device already described. Clinical trials have shown that GyneFix® mini does not increase menstrual blood loss after the first few months following insertion. Although abnormal bleeding still may occur, the low removal rate for bleeding/pain is probably due to the low incidence of concomitant complaints of pain, rendering abnormal bleeding, mostly spotting, more acceptable to the woman.

Recommendation: In case of spotting or heavy bleeding occurs the use of one of the following medications may be useful: oestrogens, oestro-progestogen, progestogen alone, NSAIDs or anti-fibrinolytic agents may be useful.

3. Expulsion of GyneFix®

With traditional IUDs, young and nulligravida/nulliparous women are particularly prone to downward migration and expulsion of the IUD. This is not the case with the GyneFix®. Long-term multicenter clinical trials with the GyneFix® showed lower expulsion rates both in parous and nulliparous women, as compared with expulsion rates of conventional copper IUDs.

Duration of action

GyneFix® is recommended for 5 years of use, or longer if recommended by the physician (women over 40).

Safety of the anchoring system

Histology studies conducted in hysterectomy specimen up to 4 years after insertion of the GyneFix® have shown that the myometrial tissue reaction at the site of the anchor is minimal (less than 1.0 mm), or even absent, supporting both the safety of the material and the safety of the implant system. No transplanted endometrial tissue could be observed within the adjacent endometrium in any of these cases studied.

Reversibility and Return to Fertility

In depth studies have shown that there is no difference in conception rates between a group who had previously used a copper IUD and a group who had not used an IUD. Fertility studies after GyneFix® removal have shown that reversal of fertility is unimpaired.

VII. Contraindications and caution (for IUDs in general)

Contraindications: pregnancy, hypoplastic uterus, abnormalities of the uterus, current pelvic inflammatory disease (PID), sexual transmitted disease (STD), postpartum endometritis or infected abortion in the past 3 months, endometrial or cervical malignancy, unexplained genital bleeding, acute cervicitis, bleeding disorders (e.g., von Willebrand's disease), genital actinomycosis, known or suspected allergy to copper, diagnosed Wilson's disease, leukaemia, chronic corticosteroid or immunosuppressive therapy.

Caution: anaemia (depending on the cause), previous ectopic pregnancy, prior uterine or tubal surgery, valvular heart disease (with antibiotic cover), heavy menstrual flow, severe menstrual cramps, fibroid uterus.

VIII. Warnings (for IUDs in general)

A. Pregnancy

1. Septic abortion

Reports on traditional IUDs have indicated an increased incidence of septic abortions, associated in some instances with septicemia, septic shock and death in patients conceiving with an IUD in place. Most of these reports have been associated with mid-trimester of pregnancy. In some cases the initial symptoms have been insidious and not easily recognized.

If a pregnancy should occur with GyneFix® in situ, the device should be removed if the string is visible even though removal may increase the chances of miscarriage. If the decision is made not to remove the device or removal proves to be difficult, or if the woman elects to continue the pregnancy, she should be warned that there may be an increased risk of abortion or preterm labour and/or sepsis and she should be followed closely. It must be noted, however, that nothing of this kind has been reported with GyneFix® so far.

2. Ectopic pregnancy

A woman who conceives while using an IUD should be carefully evaluated for the possibility of an ectopic pregnancy. Special attention should be directed to patients with delayed menses, heavy and prolonged bleeding and/or pelvic pain.

B. Pelvic Infection

An increased risk of PID associated with the use of IUDs has been reported. The risk appears to be greatest for nulliparous women who have a multiplicity of sexual partners. Studies conducted with GyneFix® in nulligravida and nulliparous women, however, have not shown any increased risk of PID in properly selected nulligravida and nulliparous women. These studies confirm that the non-prescription of IUDs to nulliparous women because of the fear of pelvic inflammatory disease (PID) and the potential for resulting infertility, (a recommendation expressed throughout the 1980s), is no longer justified. These results support the current view that it is not the age or parity of the woman but her life-style which may enhance the risk of PID. If salpingitis should occur, however, the infection could result in tubal damage and occlusion, thereby threatening future fertility. Therefore, it is recommended that women be advised to look for and report symptoms of pelvic inflammatory disease. The symptoms of pelvic infection include: development of menstrual anomaly (prolonged or heavy bleeding), abnormal vaginal discharge, abdominal or pelvic pain, dyspareunia, fever. The signs and symptoms are especially significant if they occur following the first few cycles after insertion. Appropriate aerobic and anaerobic bacteriological studies should be done and antibiotic therapy initiated promptly. If the infection does not show marked clinical improvement within 24 to 48 hours, GyneFix® should be removed and treatment reassessed on the basis of the results of culture and sensitivity tests.

C. Perforation

Partial or total perforation of the uterine wall or cervix has occurred with the use of conventional IUDs, most commonly during IUD insertion. Although perforations with GyneFix® are rare (approximately 1-2/1000), the possibility of perforation must be kept in mind during insertion and at the time of any subsequent examination. If perforation occurs, the device should be removed. Laparoscopy or laparotomy may be indicated to remove the implant. Abdominal adhesions, intestinal penetration, obstruction, and local inflammatory reaction with abscess formation and erosion of adjacent viscera may result if the implant is left in the abdominal cavity. These complications are extremely rare.

D. Effects of Copper

Additional amounts of copper available to the body from a copper IUD may precipitate symptoms in women with undiagnosed Wilson's disease. The estimated incidence of Wilson's disease is 1 in 200,000.

IX. How Supplied and How Stored

Each GyneFix® implant is sterile and delivered in a blister-pack, together with an inserter unit and uterine sound. Store dry and at room temperature.

X. Manufacturer

Contrel N.V. Incubation & Innovation Center, Technology Park, University Ghent, 9052 Zwijnaarde (Ghent), Belgium.

INSERTION INSTRUCTIONS

(See also Video INSERTION INSTRUCTIONS on www.wildemeersch.com/video/mark1.html)

Since GyneFix® represents a unique design in intrauterine contraception, physicians are cautioned that it is imperative for them to make themselves thoroughly familiar with the instructions for insertion before attempting placement of GyneFix®. Prior to insertion, the physician must review all aspects associated with the use of the GyneFix®. The patient should be given the opportunity to discuss fully any questions she may have.

A relevant medical history should be obtained to determine conditions that might influence the selection of the IUD as a method of contraception. Physical examination should include a pelvic examination and, if indicated, a Pap smear, and appropriate tests for other forms of genital disease. Pregnancy should be ruled out prior to insertion.

If the woman is anxious, the use of local intra-cervical or local/regional anaesthesia should be considered. The use of misoprostol prior to fitting of the IUD may sometimes also be useful to dilate the cervical canal. IUD providers should be aware that the popularity of the IUD could be much improved if attention is given to pain relief during IUD insertion.

After disinfection and sounding of the uterus, GyneFix® can be inserted.

The flange on the outer tube is adjusted to the measured uterine length (plus 0.5 cm to allow for the flexibility of the fundal tissue). The insertion tube is provided with centimetre marks facilitating correct placement of the flange. Note: It is not necessary to adjust the flange on the insertion tube if the provider is sufficiently familiar with the insertion technique.

Ensure alignment of the cervical canal and uterine cavity as the inserter is not very flexible. The tenaculum or atraumatic forceps or similar forceps should not be placed too high on the cervix to allow proper traction to straighten the uterine axis, especially in case of marked retroversion or anteversion.

Figure 1: Sound the uterus to obtain information on the direction and depth of the uterine cavity.

Figure 2: Adjust the flange on the insertion tube (plus 0.5 cm to allow for the flexibility of the fundal tissue).

Figure 3: Introduce the inserter until you feel that the tip making contact with the fundus. Keep the inserter in contact with the fundus by pressing the plunger forward with the thumb.

Figure 4: Move the plunger forward gently and slowly, and feel the anchor penetrating the uterine muscle.

Figure 5: Release the thread while keeping the inserter firmly against the fundus with the thumb

You can also cut the thread with a scalpel.

Figure 6: Keep traction on the uterus while removing the plunger.

Figure 7: Remove the tube carefully (removal while rotating the tube is preferred to avoid traction on the tail).

Figure 8: Trim the thread with sharp scissors. The thread can be cut in the cervix to avoid penile irritation. Alternatively the end of the thread can be pushed in the cervical canal with a long forceps, thereby making a loop. The proper positioning of the GyneFix® can be checked by measuring the distance between the upper border of the first copper tube and the serosal surface of the uterus (SS distance which is on average between 12 and 16 mm).

NOTE:

It is important not pull at the tail immediately after the procedure.

The woman should be advised to abstain from intercourse or from the use of tampons for approximately 5 days.

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